


# CPD(A) & SAGM FILLING LINE KIS



## High precision filling of blood bags

- High accuracy and filling speed of pre-assembled blood bag sets
- High filling accuracy due to mass flow metering systems.
- Automatic CIP/SIP function – ease of operation/maintenance
- Additional processing options: Seal-welding of the filling tubes, sealing the filling tube with the donor needle, vacuuming of bags, injection of provided sterile air



### Technical data



Filling volume

5 - 200 ml

Filling tolerance

100 ml +/- 1,5 ml  
50 ml +/- 1,0 ml  
20 ml +/- 0,5 ml

Output approx.  
(depending on an machine configuration)

up to 2,000 sets/h  
when both CPD(A) and SAGM filled  
  
up to 4,000 bags/h  
when only one type of fluid filled